DIGITAL MICROSCOPE	CONTENTS	PREFACE	PRODUCT DESCRIPTION	MODEL INTRODUCTION	- ES1 -	WIRELESS CAMERA PAGE AS SHOWN BELOW	OPENING THE SOFTWARE REQUIREMENT FOR "PHOTO BOOTH" SOFTWARE	CHARGING INSTRUCTIONS	TECHNICAL PARAMETERS	TROUBLESHOOTING
		English ..... 01-11		Wireless Camera: This product can be used to observe specimens in the observation mode. Long press in observation mode to enter "camera" mode. Short press "Start Recording" and short press "End Recording" again.	Windows Camera: This product can be used to observe specimens in the observation mode. Long press in observation mode to enter "camera" mode. Short press "Zoom In" or "Zoom Out" in observation mode. Long press "Image rotation" in observation mode. In other modes, move the left or right.	Windows Camera page, as shown below	Opening the software requirement for "Photo Booth" software that comes with your Apple computer.	1. Firstly, in the application directory of the Finder window, find an application named "Photo Booth".	Screen parameters: 4.3" IPS Pixels: 2.0 megapixels Magnification: 1000X Image format: JPG Video format: AVI Light source: LED Battery capacity: 1200mAh Operating time: 1.5 hours Charging time: 2 hours Supported systems: Windows vista/97/8/10/11, Mac OS X 11.0 or higher	G1: The device cannot be charged. A1: 1. Check if the power supply device is connected. 2. Turn on the main power of the power supply device (e.g., power strip, computer power).
		Deutsch ..... 12-22		PC Camera mode: Connect to a computer for TF card reading, writing, and image processing.	PC Camera mode: Connect to a computer for TF card reading, writing, and image processing.	PC Camera page, as shown below	Opening the software requirement for "Photo Booth" software that comes with your Apple computer.	2. Charging: After the charging is complete, please turn off the microscope panel. When charging is complete, the battery indicator is as shown below.	Power consumption: 2.0 megapixels Image format: JPG Video format: AVI Light source: LED Battery capacity: 1200mAh Operating time: 1.5 hours Charging time: 2 hours Supported systems: Windows vista/97/8/10/11, Mac OS X 11.0 or higher	G2: The device is frozen. A2: 1. Turn off the main power of the main unit with a three needs to force shutdown the device.
		Italiano ..... 23-33		PC Camera mode: Connect to a computer for TF card reading, writing, and image processing.	PC Camera mode: Connect to a computer for TF card reading, writing, and image processing.	PC Camera page, as shown below	Opening the software requirement for "Photo Booth" software that comes with your Apple computer.	3. After fully charged, the microscope can be used continuously for up to 1.5 hours.	PC Camera mode: Connect to a computer for TF card reading, writing, and image processing.	G3: Cannot connect to the computer. A3: Ensure the data cable is plugged in securely, replace the data cable, or restart the computer.
		Français ..... 34-44		PC Camera mode: Connect to a computer for TF card reading, writing, and image processing.	PC Camera mode: Connect to a computer for TF card reading, writing, and image processing.	PC Camera page, as shown below	Opening the software requirement for "Photo Booth" software that comes with your Apple computer.	4. After the device is completely dead, it needs to be charged with a 5V/1A charger for at least 2 hours.	PC Camera mode: Connect to a computer for TF card reading, writing, and image processing.	G4: Unable to save photos. A4: 1. If no storage card is inserted, replace it with a personal storage card. 2. If no memory card is inserted, replace it with a personal storage card. 3. Format the memory card (here: Formating all data on the card).
		Español ..... 45-55		PC Camera mode: Connect to a computer for TF card reading, writing, and image processing.	PC Camera mode: Connect to a computer for TF card reading, writing, and image processing.	PC Camera page, as shown below	Opening the software requirement for "Photo Booth" software that comes with your Apple computer.	5. After the camera is completely dead, it needs to be charged with a 5V/1A charger for at least 2 hours.	PC Camera mode: Connect to a computer for TF card reading, writing, and image processing.	G5: Image appears blurry. A5: 1. You need to rotate the focus ring to get a clearer magnified image of the object.
		日本語 ..... 56-66		PC Camera mode: Connect to a computer for TF card reading, writing, and image processing.	PC Camera mode: Connect to a computer for TF card reading, writing, and image processing.	PC Camera page, as shown below	Opening the software requirement for "Photo Booth" software that comes with your Apple computer.	6. After the device is completely dead, it needs to be charged with a 5V/1A charger for at least 2 hours.	PC Camera mode: Connect to a computer for TF card reading, writing, and image processing.	G6: Unable to save photos.
		中文 ..... 67-77		PC Camera mode: Connect to a computer for TF card reading, writing, and image processing.	PC Camera mode: Connect to a computer for TF card reading, writing, and image processing.	PC Camera page, as shown below	Opening the software requirement for "Photo Booth" software that comes with your Apple computer.	7. After the device is completely dead, it needs to be charged with a 5V/1A charger for at least 2 hours.	PC Camera mode: Connect to a computer for TF card reading, writing, and image processing.	G7: Image appears blurry.

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